



## CHIRALPLATE <sup>G</sup> special layer enantiomer separation

### Technical characteristics

- Reversed phase nano silica impregnated with Cu<sup>2+</sup> ions and a chiral selector (proline derivative)
- Separation based on ligand exchange, i.e. formation of ternary mixed-ligand complexes with the Cu(II) ions, differences in the stability of the diastereomeric complexes cause chromatographic separation

### Recommended application

- Enantiomer separation of amino acids, *N*-methylamino acids, *N*-formylamino acids,  $\alpha$ -alkylamino acids, thiazolidine derivatives, dipeptides, lactones,  $\alpha$ -hydroxycarboxylic acids


### Enantiomer separation of amino acids

MN Appl. No. 400520

Quantitative determination (remission location curves) of TLC-separated enantiomers of *tert.*-leucine:

Layer: CHIRALPLATE  
 Eluent: methanol – water (10:80, v/v)  
 Detection: dip in 0.3 % ninhydrin solution  
 quantification with scanner, 520 nm

- a) L-*tert.*-leucine  
 b) L-*tert.*-leucine + 0.1 % D-*tert.*-leucine  
 c) L-*tert.*-leucine + 1 % D-*tert.*-leucine  
 d) external reference sample



### Ordering information

Plate size [cm]	5 x 20	10 x 10	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
<b>Glass plates</b>						
Pack of [plates]	4					
CHIRALPLATE			811056		0.25 mm	UV <sub>254</sub>
Pack of [plates]	50	25	25	25		
CHIRALPLATE	811057	811059	811055	811058	0.25 mm	UV <sub>254</sub>


## SIL N-HR <sup>P</sup> unmodified standard silica layers

### Technical characteristics

- High purity silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m<sup>2</sup>/g, specific pore volume 0.75 mL/g, particle size 5–17 μm, different binder system compared to SIL G results in different separation characteristics
- A special feature of the POLYGRAM<sup>®</sup> SIL N-HR is a higher gypsum content

### Ordering information

Plate size [cm]	5 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	25		
<b>POLYGRAM<sup>®</sup> polyester sheets</b>				
SIL N-HR/UV <sub>254</sub>	804022	804023	0.20 mm	UV <sub>254</sub>



## SIL G-25 HR <sup>G</sup> special layer for aflatoxin separation

### Technical characteristics

- High purity silica 60 with gypsum and a very small quantity of a polymeric organic binder; softer than the standard silica layer, i.e. spots can be scratched and the layer absorbs faster

### Recommended application

- Aflatoxins

### Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

### Glass plates

SIL G-25 HR	809033	0.25 mm	–
SIL G-25 HR/UV <sub>254</sub>	809043	0.25 mm	UV <sub>254</sub>

## SIL G-25 Tenside <sup>G</sup> special layer for separation of surfactants

### Technical characteristics

- Silica G impregnated with ammonium sulfate

### Recommended application

- Detergents, alkanesulfonates, polyglycols

### Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

### Glass plates

SIL G-25 Tenside	810063	0.25 mm	–
------------------	--------	---------	---

## Nano-SIL PAH <sup>G</sup> special HPTLC silica layer for PAH analysis

### Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m<sup>2</sup>/g, specific pore volume 0.75 mL/g, particle size 2–10 µm
- Impregnated with caffeine, an electron acceptor for PAH analysis based on charge-transfer complexes

### Recommended application

- 6 PAHs according to German drinking water specifications (TVO) in accordance with German standard DIN 38407 part 7

### Ordering information

Plate size [cm]	10 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50		

### Glass plates

Nano-SIL PAH	811051	0.20 mm	–
--------------	--------	---------	---

Further application examples can be found online in our application database at [www.mn-net.com/apps](http://www.mn-net.com/apps)



# Layers for special TLC separations



## IONEX <sup>P</sup> special mixed layers of silica with ion exchange resins

### IONEX-25 SA-Na:

- Mixture of silica and a strongly acidic cation exchanger coated to polyester sheets

### IONEX-25 SB-AC:

- Mixture of silica and a strongly basic anion exchanger coated to polyester sheets
- Both layers contain an inert organic binder

### ✓ Recommended application

- Amino acids, e.g., in protein and peptide hydrolyzates, in seeds and fodder, in biological fluids; for racemate separation in peptide syntheses, for the separation of nucleic acid hydrolyzates, aminosugars, amino acids, antibiotics, inorganic phosphates, cations and other compounds with ionic groups

### Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

### POLYGRAM<sup>®</sup> polyester sheets

IONEX-25 SA-Na	strongly acidic cation exchanger	806013	0.20 mm	–
IONEX-25 SB-AC	strongly basic anion exchanger	806023	0.20 mm	–

## Mixed layers for TLC <sup>G</sup>

### Alox/CEL-AC-Mix-25:

- Mixed layer of aluminum oxide G and acetylated cellulose, recommended for separation of PAH

### SILCEL-Mix-25:

- Mixed layer of cellulose and silica, recommended for separation of preservatives and other antimicrobial compounds

### Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

### Glass plates

Alox/CEL-AC-Mix-25	810053	0.25 mm	–
SILCEL-Mix-25 UV <sub>254</sub>	810043	0.25 mm	UV <sub>254</sub>

Further application examples can be found online in our application database at [www.mn-net.com/apps](http://www.mn-net.com/apps)