

## Headspace screw neck vials N 18

Precision thread, not standardized



702866	702826	702826.2
10 mL	20 mL	20 mL
22.5 x 46 mm	22.5 x 75.5 mm	22.5 x 75.5 mm
clear, rounded bottom	clear, rounded bottom	amber, rounded bottom
100/PP box	100/PP box	100/PP box

## Ready assembled, magnetic screw closures N 18

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702827		N 18 magnetic screw cap, silver, center hole	Silicone blue transparent / PTFE white	45° shore A	1.5 mm	100/PE bag
702055		N 18 magnetic screw cap, silver, center hole	Silicone white / PTFE blue	45° shore A	1.5 mm	100/PE bag
702140		as above, but with closed top				
702136		N 18 magnetic screw cap, silver, center hole	Silicone white / PTFE blue, slit	45° shore A	1.5 mm	100/PE bag
702137		N 18 magnetic screw cap, silver, center hole	Butyl red / PTFE gray	50° shore A	1.5 mm	100/PE bag
702139		as above, but with closed top				
702072		N 18 magnetic screw cap, silver, center hole	Red Rubber / TEF colorless	45° shore A	1.5 mm	100/PE bag
702110		- (septum only)	Silicone white / PTFE blue	45° shore A	1.5 mm	100/PE bag
702981		- (septum only)	Silicone blue transparent / PTFE white	45° shore A	1.5 mm	100/PE bag
702138		- (septum only)	Butyl red / PTFE gray	50° shore A	1.5 mm	100/PE bag