

Screw neck vials N 24 (EPA)

For TOC/VOC analysis and sample storage



702021 / 702022	702132 / 702133	702023 / 702024	702074 / 702131
20 mL	30 mL	40 mL	60 mL
27.5 x 57 mm	27.5 x 72.5 mm	27.5 x 95 mm	27.5 x 140 mm
clear / amber, flat bottom	clear / amber, flat bottom	clear / amber, flat bottom	clear / amber, flat bottom
100 / PP box	100 / PP box	100 / PP box	100 / PP box

Container for screw neck vials N 24

REF	Description	Pack of
702517	16 position container blue, with removable divider for screw neck vials N 24 (20 mL, 30 mL, 40 mL) 130 x 130 x 102 mm, with transparent lid (suitable for freezers)	1



N 24 Screw neck vials and caps N 24

Bonded and ready assembled screw closures + plain screw caps N 24, single septa N 22

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702058		N 24 PP bonded screw cap, white, center hole	Silicone white / PTFE beige	45° shore A	3.2 mm	100 / PE bag
702059		N 24 PP bonded screw cap, white, closed top	Silicone white / PTFE beige	45° shore A	3.2 mm	100 / PE bag
702073		N 24 PP bonded screw cap, white, center hole	Red Rubber / TEF colorless	45° shore A	2.5 mm	100 / PE bag
702168		N 24 PP bonded screw cap, white, center hole, with removeable protection lid	Silicone natural / PTFE colorless	45° shore A	3.2 mm	100 / PE bag
702130		N 24 PP screw cap, white, center hole	Butyl red / PTFE gray	50° shore A	2.4 mm	100 / PE bag
702102		N 24 PP screw cap, white, closed top	Butyl red / PTFE gray	50° shore A	2.4 mm	100 / PE bag
702171		N 24 PP screw cap, white, center hole	Silicone white / PTFE beige	45° shore A	3.2 mm	100 / PE bag
702060		N 24 PP screw cap, white, center hole	no liner	-	-	100 / PE bag
702061		N 24 PP screw cap, white, closed top	no liner	-	-	100 / PE bag
702062		- (septa only)	Silicone natural / PTFE colorless	45° shore A	3.2 mm	100 / PE bag
702791		- (septa only)	Butyl red / PTFE gray	50° shore A	2.4 mm	100 / PE bag

Pre-sealed vial-closure combinations with screw neck N 24

Vial description	Closure description	Pack of	REF
Pre-sealed vials 702021: 20 mL screw neck vial N 24, 27.5 x 57 mm, clear, flat bottom	pre-screwed with 702073: N 24 PP screw cap (bonded), white, center hole, Red Rubber / TEF colorless, 45° shore A, 2.5 mm	100	702865
Pre-sealed vials 702021: 20 mL screw neck vial N 24, 27.5 x 57 mm, clear, flat bottom	pre-screwed with 702059: N 24 PP screw cap (bonded), white, closed top, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702884
Pre-sealed vials 702132: 30 mL screw neck vial N 24, 27.5 x 72.5 mm, clear, flat bottom	pre-screwed with 702058: N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702831
Pre-sealed vials 702023: 40 mL screw neck vial N 24, 27.5 x 95 mm, clear, flat bottom	pre-screwed with 702058: N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702877
Pre-sealed vials 702131: 60 mL screw neck vial N 24, 27.5 x 140 mm, amber, flat bottom	pre-screwed with 702058: N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702886
Other pre-sealed vial-closure combinations on request			

